


# Stool Testing: Microbiology Specimen Collection Guide



	Routine Stool Cultures, Giardia and Cryptosporidium EIA negative? Other Risk Factors?				
Risk	Foodborne exposure with negative cultures	Immunocompromised		Viral risk group	
Diagnostic Test	Cyclospora and Isospora Exam (CYISO) x 1-3 samples on separate days	Cyclospora and Isospora Exam (CYISO) x 1-3 samples on separate days	Microsporidia Detection (MSPORI) x1-3 samples on separate days		Immunocompromised: various agents, including Adenovirus 40/41 antigen
Collection Device	Total-Fix	Total-Fix		Sterile Container at 2-8°C	

Call Geisinger Medical  
Laboratories Client Services  
Department:  
1-800-695-6491  
with questions.

**Geisinger**

For additional information, refer to GML Test Catalog  
Click on address below  
<http://www.geisingermedicalabs.com/catalog/index.cfm>

Revised 01/14/2019  
GML stool collection.vsd