



**\*NOTE:** Gen-Probe DNA Amplified Molecular test for *C. trachomatis*/*N. gonorrhoeae* is not FDA cleared for testing on prepubescent patients ( $\leq 14$  years) or on known cases of sexual abuse. Cultures for *C. trachomatis* and *N. gonorrhoeae* are indicated in these cases. However, if collected, laboratory will perform test with above disclaimer.

**Molecular Transport Kits for *N. gonorrhoeae* & *C. trachomatis***



- 1 Use for collection and transport of swab specimens from male urethral and female endocervical specimens.
  - Amplified *N. gonorrhoeae*
  - Amplified *C. trachomatis*

Note: Use only sample collection swab included in kit.
- 2 Use for urine transport
  - Amplified *N. gonorrhoeae*
  - Amplified *C. trachomatis*
- 3 Use for collection of female vaginal specimens only.
  - Amplified *N. gonorrhoeae*
  - Amplified *C. trachomatis*

Note: Use only sample collection swab included in kit.

***Chlamydia trachomatis* molecular amplified probe:**  
<http://www.geisingermedicallabs.com/catalog/details.cfm?tid=1284>

***Chlamydia trachomatis* culture:**  
<http://www.geisingermedicallabs.com/catalog/details.cfm?tid=654>

**Collection instructions can be found at:**  
[http://www.geisingermedicallabs.com/MicroApp/molecular\\_id.shtml#chlamydia](http://www.geisingermedicallabs.com/MicroApp/molecular_id.shtml#chlamydia)

**NOTE: *Trichomonas vaginalis* testing guide information can be found at:**  
[http://www.geisingermedicallabs.com/10xEssentials/GML\\_Tvag\\_Collection.pdf](http://www.geisingermedicallabs.com/10xEssentials/GML_Tvag_Collection.pdf)

***Neisseria gonorrhoeae* molecular amplified probe:**  
<http://www.geisingermedicallabs.com/catalog/details.cfm?tid=748>

***Neisseria gonorrhoeae* culture:**  
<http://www.geisingermedicallabs.com/catalog/details.cfm?tid=1092>