

HEMOSURE[®] iFOB TEST

One-step Immunological Fecal Occult Blood Test

Patient Instructions

For assistance with these instructions contact:

- Your physician's office
 - or -
 - Hemosure technical support: 1-888-Hemosure (436-6787)
- Monday - Friday, 8:00 am to 5:30 pm, P.S.T.

STEP 1


Sample Deposit

Fill out personal information on attached label on tubes.



Attention!

- A** Lift toilet seat and position sample collection paper across the rim of the toilet bowl. Secure adhesive tabs to the sides of toilet rim. Lower the toilet seat.



- B** Make bowel movement onto collection paper.



STEP 2

Sample Collection

- A** Unscrew the purple cap from the sample collection tube. **DO NOT POUR OUT THE LIQUID.**



- B** Poke spiral applicator into stool at **6** different sites. Use only enough fecal material to cover the tip of the applicator. **DO NOT CLUMP, SCOOP, OR FILL THE TUBE.**



- C** Screw the applicator back into the tube and secure tightly.



STEP 3

Sample Return

- A** Complete the address return envelope to your doctor or laboratory.



- B** Insert sample collection tube into specimen pouch and seal. Insert specimen pouch into return envelope and seal.



- C** Return the sample pack immediately to the clinic laboratory by mail or in person.

